

LA-UR-02-2589 C.1

Approved for public release;
distribution is unlimited.

Title: Stability of Pu (VI) in synthetic WIPP brines,
chloride brines and HCl, HNO₃, HClO₄ solutions

Author(s): Patricia Paviet-Hartmann, Ningping Lu

Submitted to: DOE/Carlsbad Field Office, Carlsbad, NM



Los Alamos NATIONAL LABORATORY

Los Alamos National Laboratory, an affirmative action/equal opportunity employer, is operated by the University of California for the U.S. Department of Energy under contract W-7405-ENG-36. By acceptance of this article, the publisher recognizes that the U.S. Government retains a nonexclusive, royalty-free license to publish or reproduce the published form of this contribution, or to allow others to do so, for U.S. Government purposes. Los Alamos National Laboratory requests that the publisher identify this contribution as work performed under the auspices of the U.S. Department of Energy. Los Alamos National Laboratory strongly supports academic freedom and a researcher's right to publish; as an institution, however, the Laboratory does not endorse the viewpoint of a publication or guarantee its technical correctness.

TP 02-14-02
Revision 0
Page 1 of 1

**Los Alamos National Laboratory
EES-12
Test Plan, TP 021402, Rev.0
Stability of Pu (VI) in synthetic WIPP brines,
chloride brines and HCl, HNO₃, HClO₄ solutions**

Rev. 0
Effective Date:

Prepared by:

Patricia Paviet-Hartmann,
LANL/E-CO
Los Alamos National Laboratory
Carlsbad, NM 88220

Ningping Lu,
LANL/E-ET
Los Alamos National Laboratory
Carlsbad, NM 88220

Test Plan for the Stability of Pu(VI) in synthetic WIPP brines, chloride brines,
HCl, HNO₃, HClO₄ solutions, TP 021402, Rev. 0

TP 02-14-02
Revision 0
Page 2 of 2
1.0 APPROVAL PAGE

Author:

Patricia Paviet-Hartmann, EES-12,
LANL-PI

Date: 02-14-02

Author:

Ningping Lu, C-INC,
LANL-PI

Date: 02-14-02

LANL Technical Reviewer:

Thomas Hartmann, EES-12

Date: 02-16-02

LANL QA Reviewer:

Betty Strietelmeier, C-INC

Date:

TP 02-14-02

Revision 0

Page 3 of 3

2.0 TABLE OF CONTENTS

1.0	APPROVAL PAGE.	2
2.0	TABLE OF CONTENTS.	3
3.0	REVISION HISTORY.	4
4.0	DEFINITION OF ACRONYMS.	4
5.0	PURPOSE AND SCOPE.	7
6.0	EXPERIMENTAL PROCESS DESCRIPTION.	9
6.1	Overall Strategy and Process.	9
6.2	Planning and Quality Assurance.	19
6.3	Sample Control.	20
6.4	Data Quality Control.	21
6.5	Data Identification and Use.	21
6.6	Records, Reports, and Audits.	21
7.0	TRAINING.	22
8.0	HEALTH AND SAFETY.	22
9.0	LICENSING/PERMITTING.	23
10.0	DELIVERABLES.	23
11.0	REFERENCES.	24

3.0 REVISION HISTORY

This is the first version of this test plan (Revision 0, or Rev. 0). If it becomes necessary to revise this test plan, it will be done in accordance with the Los Alamos National Laboratory (LANL) Actinide chemistry program and it will be re-issued in accordance with LANL-EES-12 "Document Control Process." The history of any and all revisions will be described in this subsection.

4.0 DEFINITION OF ACRONYMS

CBFO: (U.S. DOE) Carlsbad Field Office.

CCA: (WIPP) Compliance Certification Application.

CEMRC: Carlsbad Environmental Monitoring Research Center

C-INC: (LANL) Chemistry Division (Division)-Isotopes and Nuclear Chemistry

DOE: (U.S.) Department of Energy.

DRC: document review and comment.

EES-12: (LANL) Environmental and Environmental Sciences Division-Carlsbad Operations.

ERDA: (U.S.) Energy Research and Development Administration.

EXAFS: Extended X ray Absorption Fine Structure

ES&H: environment safety and health.

GWB: Groundwater B, synthetic WIPP brine.

I: Ionic Strength.

ILW: Intermediate Level Waste.

LANL: Los Alamos National Laboratory.

LSC: liquid scintillation counting.

m: molal.

M: molar.

MgO: magnesium oxide.

M&TE: measuring and test equipment.

NP: (SNL) Nuclear-Waste-Management-Program procedure.

NWMP: (SNL) Nuclear Waste Management Program.

PA: performance assessment.

PI: principal investigator.

PR: purchase requisition.

Pu: plutonium (oxidation state unspecified).

TP 02-14-02

Revision 0

Page 6 of 6

Pu(III): plutonium in the +III oxidation state.

Pu(IV): plutonium in the +IV oxidation state.

Pu(V): plutonium in the +V oxidation state.

Pu(VI): plutonium in the +VI oxidation state.

Pu(VII): plutonium in the +VII oxidation state.

QA: quality assurance.

QAP: quality assurance procedure.

Rev.: revision.

SCM: (SNL) Software Configuration Manager.

SDR: Sandia-Delegated Representative.

SNL: Sandia National Laboratories.

SP: activity- or project-specific procedures.

STTP: (WIPP) Source-Term Waste Test Program.

TA: Technical Area

TP: test plan.

TRU: transuranic waste.

UV/VIS: ultraviolet/visible (light) spectrophotometry.

WIPP: Waste Isolation Pilot Plant.

WBS: Work Breakdown Structure.

WPO: (SNL) WIPP Project Office.

XRD: X-ray diffraction.

5.0 PURPOSE AND SCOPE

Because geological salt formations are considered possible sites for radioactive waste disposal, the effect of radiolysis on high-saline brine under simulated repository conditions has been the subject of numerous investigations within the last two decades (Foerster F. 1917; Allen A. 1952; Lister M. 1956; Rabideau S.W. 1957; Anbar M. 1964; Pourbaix M. 1966; Jayson G. 1973; Thornton A. 1973; Jenks G. 1975; Gilbert G. 1977; Gray W. 1985; Kim J.I. 1985; Kornilov A. 1985; Magirius S. 1985; Buppelmann K. 1986; Kim J.I. 1987; Nagar M.S. 1987; Buppelmann K. 1988; Pashalidis I. 1993; Silber H. 1994; Fukasawa T. 1996; Kelm M. 1999; Kelm M. 2000; Tandon L. 2000; Kelm M. 2001). Also of interest is radiation-chemical behavior in brines, because rock salts in geological repositories generally contain a small amount of brine inclusions that are not homogeneously distributed. Over time, these inclusions may migrate towards the waste container, which acts as a heat source. Brines will also form in geological repositories if the groundwater reaches the salt deposits (Tandon L. 2000).

In the near-field chemistry of a salt repository, the radiolytically-induced redox reactions in concentrated saline solution are of particular importance because the radiolysis of saline solutions results in the creation of oxidizing chlorine species (Anbar M. 1964; Jayson G. 1973; Thornton A. 1973; Jenks G. 1975; Gilbert G. 1977; Zansokhova A. 1977), which may oxidize actinide species to higher oxidation states

TP 02-14-02

Revision 0

Page 8 of 8

(Magirus S. 1985; Buppelmann K. 1986) but also in reducing species such as hydrogen peroxide which may reduce actinide species to lower oxidation states. Even though the number of publications concerning the influence of radiolysis by-products on the speciation of actinides is limited, when synthetic brine applications are considered, most of the studies report that hypochlorite generated by alpha irradiation (Kelm M. 1999; Kelm M. 2000) of a concentrated NaCl solution oxidizes Am (III) to Am (V), Pu (IV) to Pu (VI), and Np (IV) to Np (V) (Magirus S. 1985; Buppelmann K. 1986; Kim J.I. 1987; Buppelmann K. 1988; Pashalidis I. 1993; Silber H. 1994; Fukasawa T. 1996). Morgenstern et al. studied the kinetics of the reduction of Pu (V) O_2^+ by H_2O_2 in 1.0 M NaCl solution with time and monitored Pu(V) concentration by ultrafiltration combined with solvent extraction using TTA (Morgenstern A. 1999). The reduction of PuO_2^+ was found to be first order with respect to hydrogen peroxide concentration and inverse first order with respect to hydrogen ion concentration.

The primary objective of these experiments is to clarify the complex chemistry of actual TRU waste in contact with WIPP brines as noted in the LANL Source-Term Waste Test Program (STTP) by examining the sub-effects influencing the overall pcH-Eh system. Radiolysis is considered to be the major driving force in oxidation of actinides. Alpha- and gamma-irradiation of chloride solutions (WIPP brine) has the potential to generate secondary oxidizing agents (e.g. hypochlorite, hydrogen peroxide), and further to convert low solubility Pu (III, IV)/Am (III) species to Pu (V, VI)/Am (V) species, which are highly soluble and have a strong tendency for migration in the environment. The temporary appearance of highly mobile plutonium (VI) and (V) species in the STTP experiments might not prove to be important to the WIPP repository concept. Nevertheless, to fully understand and to predict the long-term behavior of actinides under the complex high-saline repository conditions, additional experimental approaches are necessary. The objectives of this research are:

- a- The determination of the stability of Pu (VI)-239 in synthetic WIPP brines as well as chloride brine solutions
- b- The investigation of the effect of pH on the stability of Pu (VI) in brines
- c- The investigation of the effect of H_2O_2 (10^{-2} M) on the stability of Pu (VI) in brines at pH 6, 7, 8, and 9.

d- The determination of the kinetics of the reduction of Pu(VI) by H_2O_2 and its stability in different concentrations of (1) HCl (10^{-1} M to 5M), (2) NaCl (10^{-1} M to 5M), (3) HNO_3 (10^{-1} M to 5M), and (4) HClO_4 (10^{-1} M to 5M), solutions under ambient atmosphere at room temperature.

6.0 EXPERIMENTAL PROCESS DESCRIPTION

6.1 Overall Strategy and Process

LANL personnel will use actinide speciation experiments, and characterization techniques such as UV/Vis spectrophotometry, XRD to determine the influence of alpha-radiolysis by-products on the speciation of Pu (VI) (see Section 5.0, PURPOSE AND SCOPE, above).

LANL personnel will synthesize the brine composition of interest, according to the CEMRC procedure ((Walthall M. 2002).

6.1.1 Experiments

Task 1 Stability of Pu (VI) in synthetic WIPP brine and chloride brine solutions

The major chemical in the synthetic WIPP groundwater GWB is $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$, but it also contains Na_2SO_4 , NaCl, LiCl, KCl, NaBr, and $\text{Na}_2\text{B}_4\text{O}_7 \cdot 10\text{H}_2\text{O}$. Therefore, it is necessary to compare the stability of Pu (VI) in pure $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$ solution and synthetic GWB WIPP brine. Another synthetic WIPP groundwater brine called ERDA-6, where the major constituent is NaCl, will be used to compare the stability of Pu (VI) in ERDA-6 and in a pure 5m NaCl solution. The pH values of MgCl_2 , NaCl and the two synthetic WIPP brines, ERDA-6 and GWB, are around 7, it is important to know how pH is influencing the stability of Pu (VI) in these solutions. It is also necessary to know the effect of H_2O_2 on the stability of Pu (VI) in these 4 solutions at pH relevant for the WIPP repository (i.e. pH 6, 7, 8 and 9). These experiments will provide some useful and

defensible results on Pu (VI) stability under WIPP repository conditions.

The synthetic WIPP brines will be prepared according to the procedure developed by CEMRC (Walthall M. 2002).

The experimental design is described in Table I.

Table I. Experimental Design for studying the stability of Pu (VI) in Brines

MgCl ₂ .6H ₂ O	95% GWB synthetic brine	95% ERDA-6, no NaHCO ₃	NaCl
pH 6	pH 7	pH 7	pH 7
pH6 + H ₂ O ₂	pH 7 + H ₂ O ₂	pH7 + H ₂ O ₂	pH7 + H ₂ O ₂
pH 7	pH 8	pH 8	pH 8
pH 7 + H ₂ O ₂	pH 8 + H ₂ O ₂	pH 8 + H ₂ O ₂	pH 8 + H ₂ O ₂
		pH 9	pH 9
		pH 9 + H ₂ O ₂	pH 9 + H ₂ O ₂

The experiments will be conducted in duplicate. A total of 40 samples will be analyzed by UV/Vis spectrophotometry over time (10 minutes, 1 day, 10 days, 30 days, 3 months, 6 months, 1 year). Each sample contains 3 mL of 5×10^{-4} M of Pu (VI) –239, therefore a total of 14.34 mg of Pu-239 will be used.

The ²³⁹Pu (VI) stock solution (14.34 mg) will be made up in a final matrix of dilute HCl acid (pH~2).

In a fume hood, the stock solution of Pu (VI) will be opened. With an Eppendorf pipette, the appropriate volume of stock solution will be transferred from the stock solution container to a UV/Vis cell. The volume of Pu (VI) stock solution transferred will correspond to a final activity of 5×10^{-4} M Pu(VI). Brine solution will be added to Pu (VI) to a total volume of 2.5 mL. The pH will be adjusted to 6, 7, 8, or 9 by adding the appropriate amount of NaOH solution, the proper amount of H₂O₂ will be added to the solution. The final volume will be made up to 3mL of Pu (VI) in brine in the UV/Vis

TP 02-14-02

Revision 0

Page 11 of 11

cell. The UV/Vis cell will be sealed for an indefinite period of time (~1 year). The cells will be stored in a small controlled-atmosphere box (10 inch * 10 inch). Periodically, an absorbance UV/Vis spectrum will be taken which implies transferring the sealed UV/Vis cell containing the Pu (VI) in brine (5×10^{-4} M) from the fume hood to the spectrophotometer. The lid of the spectrophotometer will be opened, the UV/Vis cell will be placed in the spectrophotometer. A spectrum will be taken over the wavelength range between 200 nm to 1100 nm using a time period of approximately 1 minute.

When the operation is finished, the UV/Vis cell containing the Pu (VI) solution will be transferred from the spectrophotometer to its original place in the controlled-atmosphere box in the fume hood.

Notes:

Preliminary tests for adjusting the pH of brine solutions have been performed on non-active brine solutions by adding drops of NaOH of different concentrations, by Dr. Mark Walthall from the Carlsbad Environmental and Monitoring Research Center, Carlsbad.

The results obtained are as follows:

3.7M $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$ - pH 7 no precipitate

3.7M $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$ - pH 8, and 9 precipitate

95% GW-Brine - pH 7, pH 8 no precipitate

95% GW-Brine - pH 9 precipitate

95% ERDA-6 with NaHCO_3 - pH 7 no precipitate

95% ERDA-6 with NaHCO_3 - pH 8, and 9 precipitate

95% ERDA-6 without NaHCO_3 - pH 7, pH 8, pH 9 no precipitate

5M NaCl - pH 7, pH 8, pH 9 no precipitate

These results were used in the conceptual design of the actual experiments (see Table 1).

Task2: Reduction of Pu(VI) by H_2O_2 in HCl solutions of varying concentrations

TP 02-14-02
Revision 0
Page 12 of 12

The reduction of Pu(VI) by H_2O_2 will be studied in HCl solutions. Table II is describing the experimental matrix.

Table II: Experimental Design for studying the stability of Pu (VI) in HCl

Sample	HCl (M) without H_2O_2 + Pu(VI) (10^{-3} M)	HCl (M) with 10^{-2} M H_2O_2 + Pu(VI) (10^{-3} M)
1	0.1	0.1
2	0.5	0.5
32	1.0	1.0
4	2.0	2.0
5	3.0	3.0
6	4.0	4.0
7	5.0	5.0

The Pu(VI) experiment will require ~5 ml of an 0.01 M Pu(VI) in 0.01 M HCl solution.

Total 14 samples will be analyzed by UV Vis spectrophotometer, along with some standards, at the following intervals:

- immediately after addition of Pu(VI)
- after 1 day
- after 10 days
- after 30 days
- 2 months
- 6 months and
- 12 months

If precipitates are formed, the solid sample will be analyzed by diffusion reflectance, XRD and if necessary by EXAFS. During the period of the experiment, ~0.1 ml of solution sample from each tube will be taken for LSC (liquid Scintillation counter) analysis. A 7mL stock solution containing 10^{-2} M Pu(VI)-239 in 0.1M HCl (pH = 2) will be prepared at the CMR building. It will be transported from the CMR building to TA 48 facility in the appropriate container via BUS-4. The stock solution containing Pu(VI) 239 will be opened in a HEPA filtered fume hood and the necessary amount of Pu-239 will be transferred in each UV Vis cell containing the HCl solution. The pH of each sample will

TP 02-14-02

Revision 0

Page 13 of 13

be determined as well as the Eh. The final volume of each sample will be 5 mL with a concentration of Pu-239 of 10^{-3} M.

The optical spectroscopy UV/Vis/NIR spectrophotometer will be used to analyze and characterize Pu containing samples. This is the method for characterizing a sample by directing a light beam (conventional light source) into the sample and detecting the response. The various oxidation states of the lighter actinides all have distinct spectra. This makes spectrophotometry uniquely valuable for determining or confirming oxidation state distributions and for assay. Capped absorption cells are used to assure no leakage around the stopper due to capillary action. These cells are filled in a hood, exterior surfaces checked for contamination, placed in a transfer container and taken to the spectrophotometer. The cells are then placed in the cell holder and spectra determined. The cells are returned to the hood via the transfer container. They are stored in a small glovebox (22 cm * 22 cm), placed under a HEPA filtered hood for one year. Once the speciation studies of Pu by UV Vis spectrophotometry are completed, the cells are emptied in a container containing vermiculite for adsorbing the Pu solution. The container will be then disposed. A waste profile form has been set up for this experiment.

Task3: Reduction of Pu(VI) by H_2O_2 in NaCl solutions of varying concentrations

The reduction of Pu(VI) by H_2O_2 will be studied in NaCl solutions. Table III is describing the experimental matrix.

Table III: Experimental Design for studying the stability of Pu (VI) in NaCl

Sample	NaCl (M) without H ₂ O ₂ + Pu(VI) (10 ⁻³ M)	NaCl (M) with 10 ⁻² M H ₂ O ₂ + Pu(VI) (10 ⁻³ M)
1	0.1	0.1
2	0.5	0.5
32	1.0	1.0
4	2.0	2.0
5	3.0	3.0
6	4.0	4.0
7	5.0	5.0

The Pu(VI) experiments will require 5 ml of an 0.01 M Pu(VI) in nanopure DI water.

Total 14 samples will be analyzed by UV Vis spectrophotometer, along with some standards, at the following intervals:

- immediately after addition of Pu(VI)
- after 1 day
- after 10 days
- after 30 days
- 2 months
- 6 months and
- 9.0 12 months

If precipitates are formed, the solid sample will be analyzed by diffusion reflectance, XRD and if necessary by EXAFS. During the period of the experiment, ~0.1 ml of solution sample from each tube will be taken for LSC (liquid Scintillation counter) analysis. A 7mL stock solution containing 10⁻² M Pu(VI)-239 in 0.1M HCl (pH = 2) will be prepared at the CMR building. It will be transported from the CMR building to TA 48 facility in the appropriate container via BUS-4. The stock solution containing Pu(VI) 239 will be opened in a HEPA filtered fume hood and the necessary amount of Pu-239 will be transferred in each UV Vis cell containing the HCl solution. The pcH of each sample will be determined as well as the Eh. The final volume of each sample will be 5 mL with a concentration of Pu-239 of 10⁻³ M.

TP 02-14-02

Revision 0

Page 15 of 15

The optical spectroscopy UV/Vis/NIR spectrophotometer will be used to analyze and characterize Pu containing samples. This is the method for characterizing a sample by directing a light beam (conventional light source) into the sample and detecting the response. The various oxidation states of the lighter actinides all have distinct spectra. This makes spectrophotometry uniquely valuable for determining or confirming oxidation state distributions and for assay. Capped absorption cells are used to assure no leakage around the stopper due to capillary action. These cells are filled in a hood, exterior surfaces checked for contamination, placed in a transfer container and taken to the spectrophotometer. The cells are then placed in the cell holder and spectra determined. The cells are returned to the hood via the transfer container. They are stored in a small glovebox (22 cm * 22 cm), placed under a HEPA filtered hood for one year. Once the speciation studies of Pu by UV Vis spectrophotometry are completed, the cells are emptied in a container containing vermiculite for adsorbing the Pu solution. The container will be then disposed. A waste profile form has been set up for this experiment.

Task 4: Reduction of Pu(VI) by H_2O_2 in $HClO_4$ solutions of varying concentrations

The Pu(VI) experiments will require 5 ml of an 0.01 M Pu(VI) in 0.1 M $HClO_4$ solution. A total of 14 samples (see Table IV) will be analyzed by UV Vis spectrophotometer, along with some standards, at the following intervals:

- immediately after addition of Pu(VI)
- after 1 day
- after 10 days
- after 30 days
- 2 months
- 6 months and
- 12 months

Table IV: Experimental Design for studying the stability of Pu (VI) in HClO₄

Sample	HClO ₄ (M) without H ₂ O ₂ + Pu(VI) (10 ⁻³ M)	HClO ₄ (M) with 10 ⁻² M H ₂ O ₂ + Pu(VI) (10 ⁻³ M)
1	0.1	0.1
2	0.5	0.5
32	1.0	1.0
4	2.0	2.0
5	3.0	3.0
6	4.0	4.0
7	5.0	5.0

If precipitates are formed, the solid sample will be analyzed by diffusion reflectance, XRD and if necessary by EXAFS. The use of hard X-ray absorption spectroscopies are ideal for determining structural data of plutonium at the solid-solution interface, as structural features such as the oxidation state, the presence or absence of the linear trans-dioxo unit, the coordination in the equatorial sphere around Pu(III-VI), and the near neighbor heavy atoms.

During the period of the experiment, ~0.1 ml of solution sample from each tube will be taken for LSC (liquid Scintillation counter) analysis.

A 7mL stock solution containing 10⁻² M Pu(VI)-239 in 0.1M HClO₄ (pH = 2) will be prepared at the CMR building. It will be transported from the CMR building to TA 48 facility in the appropriate container via BUS-4. The stock solution containing Pu(VI) 239 will be opened in a HEPA filtered fume hood and the necessary amount of Pu-239 will be transferred in each UV Vis cell containing the HClO₄ solution. The pcH of each sample will be determined as well as the Eh. The final volume of each sample will be 5 mL with a concentration of Pu-239 of 10⁻³ M.

The optical spectroscopy UV/Vis/NIR spectrophotometer will be used to analyze and characterize Pu containing samples. This is the method for characterizing a sample by directing a light beam (conventional light source) into the sample and detecting the response. The various oxidation states of the lighter actinides all have distinct spectra.

This makes spectrophotometry uniquely valuable for determining or confirming oxidation state distributions and for assay. Capped absorption cells are used to assure no leakage around the stopper due to capillary action. These cells are filled in a hood, exterior surfaces checked for contamination, placed in a transfer container and taken to the spectrophotometer. The cells are then placed in the cell holder and spectra determined. The cells are returned to the hood via the transfer container. They are stored in a small glovebox (22 cm * 22 cm), placed under a HEPA filtered hood for one year. Once the speciation studies of Pu by UV Vis spectrophotometry are completed, the cells are emptied in a container containing vermiculite for adsorbing the Pu solution. The container will be then disposed. A waste profile form has been set up for this experiment.

Task 5: Reduction of Pu(VI) by H_2O_2 in HNO_3 solutions of varying concentrations

The Pu(VI) experiments will require 5 ml of an 0.01 M Pu(VI) in 0.1 M HNO_3 solution.

A total of 14 samples (table V) will be analyzed by UV Vis spectrophotometer, along with some standards, at the following intervals:

- immediately after addition of Pu(VI)
- after 1 day
- after 10 days
- after 30 days
- 2 months
- 6 months and
- 12 months

Table V: Experimental Design for studying the stability of Pu (VI) in HNO_3

Sample	HNO_3 (M) without H_2O_2 + Pu(VI) (10^{-3} M)	HNO_3 (M) with 10^{-2} M H_2O_2 + Pu(VI) (10^{-3} M)
1	0.1	0.1
2	0.5	0.5
32	1.0	1.0
4	2.0	2.0
5	3.0	3.0
6	4.0	4.0
7	5.0	5.0

If precipitates are formed, the solid sample will be analyzed by diffusion reflectance.

During the period of the experiment, ~ 0.1 ml of solution sample from each tube will be taken for LSC (liquid Scintillation counter) analysis.

A 7mL stock solution containing 10^{-2} M Pu(VI)-239 in 0.1M HNO_3 (pH = 2) will be prepared at the CMR building. It will be transported from the CMR building to TA 48 facility in the appropriate container via BUS-4. The stock solution containing Pu(VI) 239 will be opened in a HEPA filtered fume hood and the necessary amount of Pu-239 will be transferred in each UV Vis cell containing the HNO_3 solution. The pH of each sample will be determined as well as the Eh. The final volume of each sample will be 5 mL with a concentration of Pu-239 of 10^{-3} M.

The optical spectroscopy UV/Vis/NIR spectrophotometer will be used to analyze and characterize Pu containing samples. This is the method for characterizing a sample by directing a light beam (conventional light source) into the sample and detecting the response. The various oxidation states of the lighter actinides all have distinct spectra. This makes spectrophotometry uniquely valuable for determining or confirming oxidation state distributions and for assay. Capped absorption cells are used to assure no leakage around the stopper due to capillary action. These cells are filled in a hood, exterior surfaces checked for contamination, placed in a transfer container and taken to the spectrophotometer. The cells are then placed in the cell holder and spectra determined. The cells are returned to the hood via the transfer container. They are stored

in a small glovebox (22 cm * 22 cm), placed under a HEPA filtered hood for one year. Once the speciation studies of Pu by UV Vis spectrophotometry are completed, the cells are emptied in a container containing vermiculite for adsorbing the Pu solution. The container will be then disposed. A waste profile form has been set up for this experiment.

6.2 Planning and Quality Assurance

LANL personnel will carry out the work described in this test plan in accordance with SNL/WIPP NPs, formerly known as QAPs or LANL Standard Operating Procedures (SOPs) or project-specific procedures. If project specific procedures are developed, LANL Document Control program will be used to prepare and control them as designated procedures.

Each LANL individual performing work under this test plan is responsible for achieving and maintaining quality. LANL line management is responsible for verifying quality. Independent LANL, SNL, and/or DOE/CBFO assessors are responsible for ensuring that quality is adequate for the intended use and properly documented.

LANL personnel will plan and document the processes used for the experiments described in Subsection 6.1, Overall Strategy and Process (see above), in accordance with the NP 20-2, "Scientific Notebooks." If project-specific procedures are developed, Known sources of error and uncertainty, including any uncertainty about the quality of input data, and the compatibility of data processing with any conceptual/mathematical models used at each applicable stage will also be documented in the scientific notebook. Documents to be maintained as QA records (e.g. scientific notebooks) that are not already specified in this test plan will be specified in the appropriate procedures.

LANL will control documents, as appropriate, in accordance with NP 6-2, "Document Control Process."

LANL will control processes such as performing analyses, checking routine mathematical calculations, or verifying spreadsheet calculations or utility codes to ensure that these processes are accurate, consistent, and reproducible.

LANL will develop, document, review, and approve any software developed or adapted for the WIPP project that performs data reduction or analysis, and/or supports

TP 02-14-02

Revision 0

Page 20 of 20

models, and whose results are used in support of the WIPP in accordance with NP 19-1, "Software." (LANL may document simple analyses according to the requirements specified in NP 9-1, "Analyses.") They will list software classified as "off-the-shelf" or "vendor" software, data-acquisition software, and other classifications in their QA records. They will also list software classification, version, and platform and approved users in their QA records. In addition, they will submit the list to the SNL/WIPP Software Configuration Manager (SCM) for inclusion on the WIPP Software Baseline Inventory List. (They can obtain the name and address of the current SCM from the Sandia-Delegated Representative, or SDR.)

LANL will identify conditions adverse to quality, determine the cause of these conditions, correct them, take actions to prevent recurrence, and document all of these steps.

LANL personnel working in collaboration with CEMRC will identify, label and control all samples in accordance with CEMRC CP-PROC-0003-000 "Logging Research Samples into the SID database". They will transfer all samples in accordance with CEMRC CP-PROC-001-000-100898 "Procedure: Sample Tracking and Chain of Custody" or CEMRC CP-PROC-005-000 "Analytical Request From Procedure".

6.3 Sample Control

The experiments described in Subsection 6.1, Overall Strategy and Process (see above), will require some sampling of solutions and solids for post-test characterization. LANL personnel will identify, label, and control all samples in accordance with NP 13-1, "Sample Control." They will transfer all samples in accordance with SP 13-1, "Chain of Custody."

P. Paviet-Hartmann and N. Lu, the LANL Principal Investigators (PIs) responsible for the work described in this test plan, will determine the disposition of samples after this work is completed.

6.4 Data Quality Control

LANL personnel will implement a calibration program for all of the measuring and test equipment (M&TE) used for the experiments described in Subsection 6.1, Overall Strategy and Process (see above), in accordance with NP 12-1, "Control of Measuring and Test Equipment." This M&TE calibration program will meet the requirements in NP 12-1 for: (1) receiving and testing M&TE; (2) technical operating procedures for M&TE; (3) the traceability of LANL standards to nationally recognized standards such as those from the National Institute of Standards and Technology; and (4) maintaining calibration records.

6.5 Data Identification and Use

LANL personnel will use scientific notebooks to record all of the data and conclusions obtained from the experiments described in Subsection 6.1, Overall Strategy and Process (see above), in accordance with NP 20-2, "Scientific Notebooks." Data transfer and reduction controls, control of erroneous or inadequate data, including identification, segregation and disposition of the data, and data conversion controls will be documented in the scientific notebook or in a project-specific procedure (SP).

LANL will check the results of calculations in accordance with NP 9-1, "Analyses." They will use one of the following methods: (1) separate independent calculations using the same or different analytical methods as the original calculations; (2) a check of each of the calculation steps in the original calculations, or a random (spot) check of the original calculations.

6.6 Records, Reports, and Audits

LANL personnel will consider all records providing evidence of quality, including but not necessarily limited to personnel qualification and training forms, lists of M&TE and software, technical procedures, laboratory notebooks, calibration records, and

reports, to be QA records. All of these records will be accurate, complete, identifiable, and legible.

7.0 TRAINING

P. Paviet-Hartmann and N. Lu, the LANL PIs responsible for the work described in this test plan, will ensure that training of LANL personnel is assigned and completed, and that the certification of personnel qualifications is documented according to the Hazard Control Plan LANL-EES 12 HCP. LANL personnel will complete all personnel qualification and training, and all of the documentation thereof, before beginning work.

8.0 HEALTH AND SAFETY

The work described in this test plan will require that LANL personnel use small quantities of Pu. Handling Pu, and other actinide elements in the laboratory, even in small quantities, necessitates specialized requirements for facilities, procedures, training, and oversight. These environment safety and health (ES&H) requirements are even more demanding than those necessary for most other chemical elements.

LANL has extensive experience handling Pu, and other actinide elements, and has developed comprehensive, detailed ES&H requirements and standard operating procedures for such work. Description of these requirements and procedures is beyond the scope of this test plan. All LANL personnel will work under Hazard Control Plans (HCPs) for the appropriate locations where work is performed at LANL sites. All personnel will be qualified and approved to work by the appropriate LANL line manager (Group Leaders of EES-12 and C-INC).

However, it is worth noting that, in the context of LANL's experience with Pu, and other actinide elements, the quantities of Pu LANL personnel will use for this work are small, and the procedures they will use are routine. Consequently, there are no special licensing or permitting requirements required (see 9.0, LICENSING/PERMITTING, below).

TP 02-14-02
Revision 0
Page 23 of 23

9.0 LICENSING/PERMITTING

There are no licensing or permitting requirements specific to the work described in this test plan (see 8.0, HEALTH AND SAFETY, above).

10.0 DELIVERABLES

The expected deliverables for the investigation on the “Stability of Pu (VI) in synthetic WIPP brines, chloride brines and HCl, HNO₃, HClO₄ solutions” are:

- Quaterly progress reports providing experimental updates
- Final report on “Stability of Pu (VI) in synthetic WIPP brines, chloride brines and HCl, HNO₃, HClO₄ solutions”
- Journal publication for Tasks 1, 2, 3, 4 if deemed worthwhile

11.0 REFERENCES

- Allen A., H. C., Ghormley J., Davis T., (1952). "Decomposition of water and aqueous solutions under mixed fast neutrons and gamma irradiation." J. Phys. Chem., **56**: 575.
- Anbar M., T. J. (1964). J. Phys. Chem., **68**: 3829.
- Buppelmann K., K. J. I., Lierse Ch., (1988). "The redox behavior of Pu in saline solutions under radiolysis effects." Radiochimica Acta **44/45**: 65.
- Buppelmann K., M. S., Lierse Ch., Kim J.I., (1986). "Radiolytic oxidation of Am (III) to Am (V) and Pu(IV) to Pu (VI) in saline solution." J. Less Common Metals **122**: 329.
- Foerster F., Dolch P. (1917). Z. Elektrochem., **23**: 137.
- Fukasawa T., Lierse. C., Kim J.I. (1996). "Radiolytic oxidation of Np in NaCl solutions." J. Nucl. Sc. Technol. **33**(6): 486.
- Gilbert G., Ingalls R., Swallow A., (1977). Radiat. Phys. Chem., **10**: 221.
- Gray W., Simonson S. (1985). "Gamma and alpha radiolysis of salt brines." Mat. Res. Soc. Symp. Proc. **44**: 623.
- Jayson G., Parsons B., Swallow A., (1973). J. Am. Chem. Soc., **69**: 1597.
- Jenks G., Sonder E., Bopp C., Walton J., (1975). J. Phys. Chem., **79**: 871.
- Kelm M., Bohnert. E. (2000). "Radiation chemical effects in the near field of a final disposal site I: Radiolytic products formed in concentrated NaCl solutions." Nuclear Technol. **129**: 119.
- Kelm M., Bohnert E. (2001). Products formed from alpha radiolysis of chloride brines. Karlsruhe, Germany.
- Kelm M., P. I., Kim J.I. (1999). "Spectroscopic investigation on the formation of hypochlorite by alpha radiolysis in concentrated NaCl solutions." Applied radiation and isotopes **51**: 637.
- Kim J.I., Apostolidis C., Buckau G., Bueppelmann K., Kanellakopulos B., Lierse Ch., Magirius S., Stumpe R., Hedler I., Rahner Ch., Stower W., (1985). Chemical behavior of Np, Pu and Am in various brine solutions. Munich, Technische Universitat Munchen.

- Kim J.I., Lierse C., Buppelmann K., Magirus S., (1987). "Radiolitically induced oxidation reactions of actinide ions in concentrated salt solutions." Mat. res. Soc. Symp. Proc. **84**: 603.
- Kornilov A., Vasilev V., Frolov A., (1985). "Investigation of the influence of alpha irradiation on the valence states of the actinides VI. Rediction of Am(VI) in perchlorate solutions." Sov. Radiochem. **27**: 48.
- Lister M., Petterson R. (1956). "Decomposition of sodium hypochlorite the uncatalyzed reaction." Can. J. Chem., **34**: 465.
- Magirus S., Carnall W., Kim J.I., (1985). "Radiolytic oxidation of Am(III) to Am(V) in NaCl solutions." Radiochimica Acta **38**: 29.
- Morgenstern A., Choppin G. (1999). "Kinetics of the reduction of Pu(V)O_2^+ by hydrogen peroxide." Radiochimica Acta **86**: 109.
- Nagar M.S., Natarajan P. R. (1987). "Studies on radiolytic oxidation of plutonium(III) in hydrochloric acid medium." Radiation Physics and Chemistry **30**: 169.
- Pashalidis I., Kim J. I., Lierse Ch., Sullivan J., (1993). "The chemistry of Pu in concentrated aqueous NaCl solution: Effects of alpha self-radiolysis and the interaction between hypochlorite and dioxoplutonium (VI)." Radiochimica Acta **60**: 99.
- Pourbaix M. (1966). Atlas of electrochemical equilibria in aqueous solution. London, Pergamon Press: 590.
- Rabideau S.W. (1957). J. Am. Chem. Soc., **79**: 6350.
- Silber H., Nitsche H., Gatti R., Gehmecker H., Feige G., Bucher J., Edelstein N., (1994). "The effects of radiolysis upon speciation and solubility of Np in brine solutions." Radiochim. Acta **66/67**: 15.
- Tandon L. (2000). Radiolysis of Salts and Long-Term Storage Issues for Both Pure and Impure PuO_2 Materials in Plutonium Storage Containers, Los Alamos National Laboratory, LA-13725-MS
- Thornton A., Lawrence G. (1973). J. Chem. Soc., Dalton Trans.,: 804.
- Walthall M. (2002). CEMRC Brine Preparation Procedure. Carlsbad, NM, CEMRC.
- Zansokhova A., Kabakchi S., Pikaev A., (1977). Radiokhimiya **11**: 67.